Insert de détartrage supragingival - Supragingival scaling tip

Supragingival Scaler Spitze - Insert de detartraje supragingival
Inserți per detartraggio sopragengivale - Υπερουσιακό ξέστρω αποτρύγγωσης

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 **Indications:** The Scorpion scaler tip is perfectly designed for teeth scaling. The use is restricted to competent and qualified practitioners (dental surgeons, hygienists). Due to an improved vibration, power on the scaler equipment can be reduced by the dentist from 40 % in comparison with a classical tip. This offers a better comfort for the patient with an unchanged efficiency for the dentist. Screw the tip on the

handpiece using the torque wrench. Leave the insert remove the tarter without applying pressure. Maintain irrigation during treatment. The visible wear indicator, a strong but thin titanium coat, can be seen at the sharp end of the tip with its golden colour. **Disinfection – Cleaning – Sterilization:** The Scorpion tip is supplied non-sterile and must be disinfected, cleaned and sterilized before each use. Only use solutions with proven efficiency for cleaning and disinfection of instruments (e.g. DGHM, FDA or CE marked). 1) Immediately after the invasive procedure, immerse the device in a decantation bath according to the manufacturer's instructions (concentration, time, temperature). An ultrasonic bath may also be used. 2) Clean the device with running water until free of visible debris and disinfectant residues. 3) Wash the device manually using a soft brush and a cleaning solution according to the manufacturer's instructions (concentration, time, temperature). 4) Rinse thoroughly with running water to remove visible dirt and cleaning product residues. 5) Dry the device with a soft lint-free tissue or with filtered compressed air. 6) Visually check for dirt (and corrosion). 7) Place the device in proper packaging for medical devices sterilization according to ISO 11607 – EN 868. 8) Sterilize the device using autoclave (steam sterilization) at the temperature of 134 °C for 18 minutes. **Note:** Do not use on implants and restorative materials. Dispose of the tip according to the applicable regulations regarding medical waste with an infectious risk. The use and the application of the product are fully achieved under the responsibility of the dental practitioner. No liability will be accepted by the manufacturer in case of damage.

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