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
Indications: The CLiP is indicated for the elimination of the plaque and tartar on dental implants and prosthetic restorations. The use is restricted to competent and qualified practitioners (dental surgeons, hygienists). It is designed to work exclusively with the Scorpion-i ultrasonic tip and permits cleaning without damaging the fragile surface during scaling.

Fixing the CLiP onto the insert: Place the tip of the scaler into one of the CLiP's in the CLiPPER (fig. 1). Press the tip firmly into the CLiP (minimum 3 kg) to ensure good fixation (fig. 2). Remove the insert with the CLiP attached and manually check that the CLiP is securely fixed.

Cleaning implants and restorations (fig. 3): Ensure the power setting on the ultrasonic generator is set from "minimum" to 20% max. (periodontology), or Stage 1 for pneumatic handpieces. Set the water spray to a medium setting. When scaling, work with aspiration to ensure that your patient does not swallow any debris and ensure that the CLiP is parallel to the surface being cleaned. Do not press the CLiP against the surface.

To remove the CLiP: Place the tip with the CLiP into the groove at the end of the CLiPPER (fig. 4). Visualize if the CLiP is worn thanks to the witness line (fig. 5 and 6). The CLiP can be removed from the tip by pulling hard against the wear guide at the side of the CLiPPER (fig. 4).

Disinfection - Cleaning - Sterilization: The CLiPs, as well as the tips, are supplied non-sterile and must be disinfected, cleaned and sterilized before each use. Only use solutions with proved efficiency for cleaning and disinfection of instruments (e.g. DGHM, FDA or CE marked). 1) Immediately after the invasive procedure, immerse in a decontamination bath according to the manufacturer's instructions (concentration, time, temperature). An ultrasonic bath may also be used. 2) Clean with running water until free of visible debris and disinfectant residues. 3) Wash manually using a soft brush and a cleaning solution according to the manufacturer's instructions (concentration, time, temperature). 4) Rinse thoroughly with running water to remove visible dirt and cleaning product residues. 5) Dry with a soft lint-free tissue or with filtered compressed air. 6) Visually check for dirt. 7) Place in proper packaging for medical devices sterilization according to ISO 11607 - EN 868. 8) Sterilize using autoclave (steam sterilization) at the temperature of 134 °C for 18 minutes.

Elimination: The life of the CLiP depends upon the time it is in use. Dispose of the CLiP according to the applicable regulations regarding medical waste with an infectious risk.

Important : The CLiP is designed to work exclusively with Scorpion i-inserts (I-SA, I-EM, I-SI, I-KP, I-KA, I-AIR, I-AIRw). Using the CLiP with any other tips will result in the CLiP falling off when in use. When fixing the CLiP to the tip only use the CLiPPER to ensure secure fixation. Do not try to fix a CLiP to the tip without the CLiPPER for example by hand. The correct use and application of the CLiP is the responsibility of the user. Scorpion does not accept any responsibility in any case of misuse.

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- Instruction for use
- Bedienungsanleitung
- Instrucciones de empleo
- Istruzioni per l'uso
- Οδηγίες Χρήσης
- Instruções de uso
- Instruktationer
- Návod na použití
- Gebruiksaanwijzing
- Инструкции за употреба
- 使用手册
- Kullanım için talimatlar

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Formal: CLIP will be formal at first place and catalog in the CLP's and prosthetic restorations. It can be removed at the following procedure: (fig. 3). Push the insert into the CLP until it is firmly fixed. Do not force the insert into the CLP. **Preparation of CLIP for use:** Insert the tip of the insert into the CLP, but do not force the insert into the CLP (minimum 3 kg). Press the insert into the CLP correctly to ensure a good fit (fig. 2). Remove the insert with the CLP attached and check manually that the CLP is firmly fixed. **Cleaning, disinfection and sterilization (fig. 3):** Adjust the power of the ultrasonic generator to 20% maximum (Parodontology mode) and position 1 of the pneumatic handle. Immerse the spray of the insert in a medical bath. Stop and start the ablation, checking to avoid the patient accidentally touching their teeth. Make sure that the CLP is parallel to the surface to be cleaned. Do not press the CLP against the surface. **Removal of the CLP:** Place the insert with the CLP on the scale at the end of the CLIPPER (fig. 4). Check the wear of the CLP using the control line (fig. 5 and 6). The CLP can be removed from the insert by pushing it firmly against the wear tip on the side of the CLIPPER (fig. 4). **Disinfection - Cleaning - Sterilization:** CLP's (like the insert) are supplied non-sterile and must be disinfected, cleaned and sterilized before each use. Use only solutions with an effective concentration for cleaning and disinfection of instruments (for example DGHM, FDA or CE marking). 1) Immediately after use, immerse the device in a disinfection bath according to the manufacturer's instructions (concentration, time, temperature). It can also be used for ultrasonic cleaning. 2) Rinse thoroughly with water until all visible debris and disinfectant residues are removed. 3) Wash by hand with a soft brush and a cleaning solution according to the manufacturer's instructions (concentration, time, temperature). 4) Rinse thoroughly with running water and dry with a soft cloth or filtered compressed air. 5) Check visually for any debris. 6) Place the device in a suitable container for the sterilization of medical devices according to ISO 11607 - EN 868. 8) Sterilize in autoclave (vapor - sterilization with color change) at a temperature of 134 ° C for 18 minutes. **Disposal:** The life of the CLP depends on use and wear (movement) during use. Dispose of the CLP in accordance with the rules for disposal of medical waste and infection risk. **Important:** The CLP is designed to work exclusively with Scorpion i inserts (SA, EM, SI, KP, KA, IAR, IARW). Use of the CLP with another type of insert may cause it to fall during use. To fix the CLP to the insert, use only the CLIPPER to ensure a correct fit. Do not try to fix the CLP to the insert without the CLIPPER, for example by hand. The user is responsible for the use and correct installation of the CLP. Scorpion declines any liability in the event of incorrect use.

Using instructions: " CLIP" is intended for the removal of plaque and tartar deposits from dental implants and prosthetic restorations. The use of CLIP is limited exclusively to qualified and experienced dental practitioners (oral surgeons, hygienists). It is used exclusively on Scorpion i ultrasonic inserts, allowing safe cleaning without the risk of damaging the surface during the ablation. **Preparation of CLIP for use:** Insert the tip of the insert into the CLP, but do not force the insert into the CLP (minimum 3 kg). Press the insert into the CLP correctly to ensure a good fit (fig. 2). Remove the insert with the CLP attached and check manually that the CLP is firmly fixed. **Cleaning, disinfection and sterilization (fig. 3):** Adjust the power of the ultrasonic generator to 20% maximum (Parodontology mode) and position 1 of the pneumatic handle. Immerse the spray of the insert in a medical bath. Stop and start the ablation, checking to avoid the patient accidentally touching their teeth. Make sure that the CLP is parallel to the surface to be cleaned. Do not press the CLP against the surface. **Removal of the CLP:** Place the insert with the CLP on the scale at the end of the CLIPPER (fig. 4). Check the wear of the CLP using the control line (fig. 5 and 6). The CLP can be removed from the insert by pushing it firmly against the wear tip on the side of the CLIPPER (fig. 4). **Disinfection - Cleaning - Sterilization:** CLP's (like the insert) are supplied non-sterile and must be disinfected, cleaned and sterilized before each use. Use only solutions with an effective concentration for cleaning and disinfection of instruments (for example DGHM, FDA or CE marking). 1) Immediately after use, immerse the device in a disinfection bath according to the manufacturer's instructions (concentration, time, temperature). It can also be used for ultrasonic cleaning. 2) Rinse thoroughly with water until all visible debris and disinfectant residues are removed. 3) Wash by hand with a soft brush and a cleaning solution according to the manufacturer's instructions (concentration, time, temperature). 4) Rinse thoroughly with running water and dry with a soft cloth or filtered compressed air. 5) Check visually for any debris. 6) Place the device in a suitable container for the sterilization of medical devices according to ISO 11607 - EN 868. 8) Sterilize in autoclave (vapor - sterilization with color change) at a temperature of 134 ° C for 18 minutes. **Disposal:** The life of the CLP depends on use and wear (movement) during use. Dispose of the CLP in accordance with the rules for disposal of medical waste and infection risk. **Important:** The CLP is designed to work exclusively with Scorpion i inserts (SA, EM, SI, KP, KA, IAR, IARW). Use of the CLP with another type of insert may cause it to fall during use. To fix the CLP to the insert, use only the CLIPPER to ensure a correct fit. Do not try to fix the CLP to the insert without the CLIPPER, for example by hand. The user is responsible for the use and correct installation of the CLP. Scorpion declines any liability in the event of incorrect use.

Using instructions: " CLIP" is intended for the removal of plaque and tartar deposits from dental implants and prosthetic restorations. The use of CLIP is limited exclusively to qualified and experienced dental practitioners (oral surgeons, hygienists). It is used exclusively on Scorpion i ultrasonic inserts, allowing safe cleaning without the risk of damaging the surface during the ablation. **Preparation of CLIP for use:** Insert the tip of the insert into the CLP, but do not force the insert into the CLP (minimum 3 kg). Press the insert into the CLP correctly to ensure a good fit (fig. 2). Remove the insert with the CLP attached and check manually that the CLP is firmly fixed. **Cleaning, disinfection and sterilization (fig. 3):** Adjust the power of the ultrasonic generator to 20% maximum (Parodontology mode) and position 1 of the pneumatic handle. Immerse the spray of the insert in a medical bath. Stop and start the ablation, checking to avoid the patient accidentally touching their teeth. Make sure that the CLP is parallel to the surface to be cleaned. Do not press the CLP against the surface. **Removal of the CLP:** Place the insert with the CLP on the scale at the end of the CLIPPER (fig. 4). Check the wear of the CLP using the control line (fig. 5 and 6). The CLP can be removed from the insert by pushing it firmly against the wear tip on the side of the CLIPPER (fig. 4). **Disinfection - Cleaning - Sterilization:** CLP's (like the insert) are supplied non-sterile and must be disinfected, cleaned and sterilized before each use. Use only solutions with an effective concentration for cleaning and disinfection of instruments (for example DGHM, FDA or CE marking). 1) Immediately after use, immerse the device in a disinfection bath according to the manufacturer's instructions (concentration, time, temperature). It can also be used for ultrasonic cleaning. 2) Rinse thoroughly with water until all visible debris and disinfectant residues are removed. 3) Wash by hand with a soft brush and a cleaning solution according to the manufacturer's instructions (concentration, time, temperature). 4) Rinse thoroughly with running water and dry with a soft cloth or filtered compressed air. 5) Check visually for any debris. 6) Place the device in a suitable container for the sterilization of medical devices according to ISO 11607 - EN 868. 8) Sterilize in autoclave (vapor - sterilization with color change) at a temperature of 134 ° C for 18 minutes. **Disposal:** The life of the CLP depends on use and wear (movement) during use. Dispose of the CLP in accordance with the rules for disposal of medical waste and infection risk. **Important:** The CLP is designed to work exclusively with Scorpion i inserts (SA, EM, SI, KP, KA, IAR, IARW). Use of the CLP with another type of insert may cause it to fall during use. To fix the CLP to the insert, use only the CLIPPER to ensure a correct fit. Do not try to fix the CLP to the insert without the CLIPPER, for example by hand. The user is responsible for the use and correct installation of the CLP. Scorpion declines any liability in the event of incorrect use.

Using instructions: " CLIP" is intended for the removal of plaque and tartar deposits from dental implants and prosthetic restorations. The use of CLIP is limited exclusively to qualified and experienced dental practitioners (oral surgeons, hygienists). It is used exclusively on Scorpion i ultrasonic inserts, allowing safe cleaning without the risk of damaging the surface during the ablation. **Preparation of CLIP for use:** Insert the tip of the insert into the CLP, but do not force the insert into the CLP (minimum 3 kg). Press the insert into the CLP correctly to ensure a good fit (fig. 2). Remove the insert with the CLP attached and check manually that the CLP is firmly fixed. **Cleaning, disinfection and sterilization (fig. 3):** Adjust the power of the ultrasonic generator to 20% maximum (Parodontology mode) and position 1 of the pneumatic handle. Immerse the spray of the insert in a medical bath. Stop and start the ablation, checking to avoid the patient accidentally touching their teeth. Make sure that the CLP is parallel to the surface to be cleaned. Do not press the CLP against the surface. **Removal of the CLP:** Place the insert with the CLP on the scale at the end of the CLIPPER (fig. 4). Check the wear of the CLP using the control line (fig. 5 and 6). The CLP can be removed from the insert by pushing it firmly against the wear tip on the side of the CLIPPER (fig. 4). **Disinfection - Cleaning - Sterilization:** CLP's (like the insert) are supplied non-sterile and must be disinfected, cleaned and sterilized before each use. Use only solutions with an effective concentration for cleaning and disinfection of instruments (for example DGHM, FDA or CE marking). 1) Immediately after use, immerse the device in a disinfection bath according to the manufacturer's instructions (concentration, time, temperature). It can also be used for ultrasonic cleaning. 2) Rinse thoroughly with water until all visible debris and disinfectant residues are removed. 3) Wash by hand with a soft brush and a cleaning solution according to the manufacturer's instructions (concentration, time, temperature). 4) Rinse thoroughly with running water and dry with a soft cloth or filtered compressed