	Insert de détartrage supragingival - Supragingival scaling tip Supragingival Scaler Spitze - Insert de detartraje supragingival Inserti per detartraggio sopragengivale - Υπερουβικό ξέστορο αποτρίγωσης 龈上工作尖 - Supragingivalnių apnašų šalinimo antgalis													
		<table><tr><td>BSA</td><td>Satelec - NSK Varios - DTE</td></tr><tr><td>BEM</td><td>EMS - Woodpecker</td></tr><tr><td>BSI</td><td>Sirona Sirosonic</td></tr><tr><td>BKP</td><td>Kavo Piezolux Sonosoft - MECTRON</td></tr><tr><td>AIR</td><td>Kavo Sonicflex 2000-2003</td></tr><tr><td>AIRw</td><td>W&H Synex-Proxoco NSK Ti-MAX S970-AS2000</td></tr></table>	BSA	Satelec - NSK Varios - DTE	BEM	EMS - Woodpecker	BSI	Sirona Sirosonic	BKP	Kavo Piezolux Sonosoft - MECTRON	AIR	Kavo Sonicflex 2000-2003	AIRw	W&H Synex-Proxoco NSK Ti-MAX S970-AS2000
BSA	Satelec - NSK Varios - DTE													
BEM	EMS - Woodpecker													
BSI	Sirona Sirosonic													
BKP	Kavo Piezolux Sonosoft - MECTRON													
AIR	Kavo Sonicflex 2000-2003													
AIRw	W&H Synex-Proxoco NSK Ti-MAX S970-AS2000													
Dispositif médical Medical device Medizinprodukt Producto sanitario Dispositivo médico Ιατρική συσκευή 醫療裝置 Medicinos priemonė	Adaptabilité Compatibilité Compatibilità Compatibilidad Compatibilità Συμβατότητα 兼容性 Suderinamumas													

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Indications: The Scorpion scaler tip is perfectly designed for teeth scaling. The use is restricted to competent and qualified practitioners (dental surgeons, hygienists). Due to an improved vibration, power on the scaler equipment can be reduced by the dentist from 40 % in comparison with a classical tip. This offers a better comfort for the patient with an unaltered efficiency for the dentist. Screw the tip on the handpiece using the torque wrench. Leave the insert remove the tartar without applying pressure. Maintain irrigation during treatment. The tip should be replaced when the gold colored Titanium Nitride usage indicator has all worn away. **Disinfection – Cleaning – Sterilization:** The Scorpion tip is supplied non-sterile and must be disinfected, cleaned and sterilized before each use. Only use solutions with proved efficiency for cleaning and disinfection of instruments (e.g. DGHM, FDA or CE marked). 1) Immediately after the invasive procedure, immerse the device in a decontamination bath according to the manufacturer's instructions (concentration, time, temperature). An ultrasonic bath may also be used. 2) Clean the device with running water until free of visible debris and disinfectant residues. 3) Wash the device manually using a soft brush and a cleaning solution according to the manufacturer's instructions (concentration, time, temperature). 4) Rinse thoroughly with running water to remove visible dirt and cleaning product residues. 5) Dry the device with a soft lint-free tissue or with filtered compressed air. 6) Visually check for dirt (and corrosion). 7) Place the device in proper packaging for medical devices sterilization according to ISO 11607 - EN 868. 8) Sterilize the device using autoclave (steam sterilization) at the temperature of 134 ° C for 3 - 18 minutes. **Note:** Do not use on implants and restorative materials. Dispose of the tip according to the applicable regulations regarding medical waste with an infectious risk. The use and the application of the product are fully achieved under the responsibility of the dental practitioner. We decline all responsibility in case of bad use of the product. Any serious incident occurring in connection with the device must be notified to the manufacturer and to the local competent authority.

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