




Universal  
scaler tip




Visual wear  
indicator

The Scorpion insert tip is coated with a 2 µm thick and strong titanium nitride layer. To maintain efficiency and for patient comfort the insert has to be changed when the titanium coating has disappeared

Very thin




Scorpion tip



Other tip

Scorpion tip reaches interproximal spaces without difficulty

Very powerful



Scorpion tip



Other tip

A strong cavitation is guaranteed even with a low power scaling device

CLiP  
& INSERT-i

Safe Implant Cleaner

Clean implants and restorations without changing your tip !



Clip x5

+




Insert-ī

+



Dynamometric key x 1

+



Clipper

## EMS - Woodpecker



## Satelec - NSK varios - DTE



## Sirona Sirosonic



## Kavo PiezoLux - SonoSoft - Mectron



## Kavo SONICflex 2000-2003



## W&H Synea-Proxeo NSK Ti-Max S970 - AS2000



## Kavo PiezoLED Kavo PIEZOsoft



## LM - Amdent



Made in France  
**SCORPION**  
 4 rue Enrico Fermi - 63540 Romagnat  
 FRANCE  
 E-mail: [contact@scorpion.fr](mailto:contact@scorpion.fr)  
 Web: [www.scorpion.fr](http://www.scorpion.fr)

