

SHRILL NOISES REDUCER

DENTAL

MAGi(C)  
SONIC





## SOUND PRESSURE DAMPER

**MAGIC SONIC** consists of 3 resonators which reduce the sound pressure in the frequency range of around 6000 Hz.

**The sound environment and the perception of speech sounds remain unchanged for the comfort of dentists, hygienists, assistants and patients.**

**MAGIC SONIC** does not create an occlusion effect (stuffy ear sensation).



PATENT

Made in France



The box contains :



**MAGIC  
SONIC**

x 4

(2 pairs)

